



Product Datasheet

PDZD8 antibody (orb1959)


Description	PDZD8 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PDZD8 (901-1000/1154 aa)
Target	PDZD8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	128 kDa
Uniprot ID	Q8NEN9
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Sample:
HeLa(Human) Cell
Lysate at 30 ug...



Paraformaldehyde-
fixed, paraffin
embedde...



Sample: Panc-
1(Human) Cell
Lysate at 30 ...