



Product Datasheet

ACTA2 antibody (orb195993)


Description	ACTA2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SM alpha-actin. Please contact us for the exact immunogen sequence. The peptide is available as orb394004.
Target	ACTA2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q9JI91 , P35609 , D3ZCV0
Entrez	88
Dilution Range	IF/ICC: 1:500, WB: 1:200, IHC-P: 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of pig larg...



IF analysis of pig large intestines tiss...



IF image of rat heart tissue using anti-...