

## Product Datasheet

### VEGF121 Protein, Human, Recombinant (His & Avi) (orb1960085)

|                        |                                                                                                                                                                                                                                |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | VEGF121 Protein, Human,<br>Recombinant (His & Avi)                                                                                                                                                                             |
| <b>Note</b>            | For research use only                                                                                                                                                                                                          |
| <b>Purity</b>          | 98.00%                                                                                                                                                                                                                         |
| <b>MW</b>              | The protein has a predicted MW of 17 kDa. Due to glycosylation, the protein migrates to 18 kDa and 22-25 kDa under reduced (R) condition, 30 kDa and 32-40 kDa under Non reducing (N) condition based on Tris-Bis PAGE result. |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                                 |