

Product Datasheet

Eotaxin 3/CCL26 Antibody (orb196237)

Description

Eotaxin 3/CCL26 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

E.coli-derived human Eotaxin 3 recombinant protein (Position: T24-L94).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for Eotaxin 3 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

10648 MW

Uniprot ID
[Q9Y258](#)
Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By HeatWestern blot, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



WB analysis of Eotaxin 3 using anti-Eotaxin 3 antibody.



IHC analysis of Eotaxin 3 using anti-Eotaxin 3 antibody.



IHC analysis of Eotaxin 3 using anti-Eotaxin 3 antibody.