




Product Datasheet

CXCR3 Antibody (orb196249)

Description		CXCR3 Antibody	
Species/Host	Rabbit		
Reactivity	Human, Rat		
Conjugation	Unconjugated		
Tested Applications	IHC, WB		
Immunogen	E.coli-derived human CXCR3 recombinant protein (Position: M1-L368). Human CXCR3 shares 86% amino acid (aa) sequence identity with both mouse and rat CXCR3.		
Form/Appearance	Lyophilized		
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.		
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.		
Note	For research use only		
Application notes	WB: The detection limit for CXCR3 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.		
Isotype	Rabbit IgG	IHC analysis of CXCR3 using anti-CXCR3 a...	
Clonality	Polyclonal		
MW	40660 MW		
Uniprot ID	P49682		
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Rat, By HeatWestern blot, 0.1-0.5µg/ml, Human		
Expiration Date	12 months from date of receipt.		