

Product Datasheet

Ubiquitin/UBB Antibody (orb196292)

Description

Ubiquitin/UBB Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human Ubiquitin recombinant protein
(Position: M77-G152).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of
500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months. Avoid
repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for Ubiquitin is approximately
0.25ng/lane under reducing conditions. Tested Species: In-
house tested species with positive results. By Heat: Boiling
the paraffin sections in 10mM citrate buffer, pH6.0, for
20mins is required for the staining of formalin/paraffin
sections. Other applications have not been tested. Optimal
dilutions should be determined by end users. . Add 0.2ml of
distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW




25762 MW

Uniprot ID
P0CG47
Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat,
Immunohistochemistry (Paraffin-embedded Section), 0.5-
1µg/ml, Human, Mouse, Rat, By Heat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Human

Expiration Date

12 months from date of receipt.

kDa 1 2 3 4 5 6
72-

WB analysis
of Ubiquitin
using anti-
Ubiquitin

IF analysis of
Ubiquitin
using anti-
Ubiquitin

IHC analysis
of Ubiquitin
using anti-
Ubiquitin