



Product Datasheet

Human IGF1R/Igf1 Receptor ELISA Kit (orb196294)

Description	Human IGF1R/Igf1 Receptor ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E31-N932
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	154793 MW
Uniprot ID	P08069
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml.
Expiration Date	Please enquire.



A typical standard curve for the Human I...



Human IGF1R PicoKine ELISA Kit standard ...