


Product Datasheet

Human IL-17RB ELISA Kit (orb196295)

Description	Human IL-17RB ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R18-K286
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	55885 MW
Uniprot ID	Q9NRM6
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human IL-17RB PicoKine ELISA Kit standard...