


Product Datasheet

Human IL-17C ELISA Kit (orb196296)

Description	Human IL-17C ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H19-V197
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	21765 MW
Uniprot ID	Q9P0M4
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 125 pg/ml.
Expiration Date	Please enquire.



Human IL-17C PicoKine ELISA Kit standard...