


Product Datasheet

Human SPARCL1/SPARC-like protein 1 ELISA Kit (orb196299)

Description	Human SPARCL1/SPARC-like protein 1 ELISA Kit
Reactivity	Human
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I17-F664
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	75208 MW
Uniprot ID	Q14515
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 100 ng/ml.
Expiration Date	Please enquire.



Human
SPARCL1
PicoKine
ELISA Kit
standar...