

Product Datasheet

SCGB2A2 antibody (orb19630)

Description

Goat polyclonal antibody to SCGB2A2

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

Mammaglobin A / SCGB2A2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Breast shows staining of cytoplasm in some acinar cells.. Recommended concentration, 4-6ug/ml.WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 10.5kDa according to NP_002402.1. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

NPQVSKTEYKEL

Clonality

Polyclonal

MW

10.5

Entrez


4250

Dilution Range

ELISA: 1:32000, IHC-P: 4-6µg/ml

Expiration Date

12 months from date of receipt.



orb19630
(4ug/ml)
staining of
paraffin e...