


Product Datasheet

Human SDC1 / Syndecan-1 / CD138 ELISA Kit (orb196303)

Description	Human SDC1 / Syndecan-1 / CD138 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q23-G254
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	32462 MW
Uniprot ID	P18827
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 68-5651 pg/ml.
Expiration Date	Please enquire.



Human SDC1/Syndecan-1 PicoKine ELISA Kit...