


Product Datasheet

Mouse Neuropilin-1 ELISA Kit (orb196306)

Description	Mouse Neuropilin-1 ELISA Kit
Reactivity	Mouse
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F22-P856
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	103000 MW
Uniprot ID	P97333
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 1000 ng/ml.
Expiration Date	Please enquire.



Mouse
Neuropilin-1
PicoKine
ELISA Kit
st...