


Product Datasheet

Rat VEGFR3/FLT4 ELISA Kit (orb196310)

Description	Rat VEGFR3/FLT4 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y25-D770
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	153317 MW
Uniprot ID	Q91ZT1
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 1 ng/ml.
Expiration Date	Please enquire.



Rat
VEGFR3/FLT4
PicoKine
ELISA Kit
stand...