


Product Datasheet

Human Livin ELISA Kit (orb196311)

Description	Human Livin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-S298
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and cell lysates.
Uniprot ID	Q96CA5
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 1.5 ng/ml.
Expiration Date	Please enquire.



Human Livin
PicoKine
ELISA Kit
standard ...