

Product Datasheet

GALNS antibody (orb19640)

Description	Goat polyclonal to GALNS
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GALNS
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 60kDa band observed in Human Bone Marrow and Placenta lysates (calculated MW of 58.0kDa according to NP_000503.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	TTHNLEDHTKLP
Clonality	Polyclonal
MW	58
Entrez	2588
Dilution Range	ELISA: 1:64000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



orb19640
(0.01ug/ml)
staining of
Human B...