


Product Datasheet

ARSA antibody (orb19647)


Description	Goat polyclonal antibody to ARSA
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	Arylsulfatase A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55kDa band observed in Mouse Testis lysates (calculated MW of 53.6kDa according to human NP_000478.2 and 53.7kDa according to Mouse NP_033843.2). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	YDLSKDPGENYN
Clonality	Polyclonal
MW	53.6; 44.9
Entrez	410
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Testis ly...



Immunohistochemical staining of Human Co...



Immunohistochemical staining of Human Te...