


Product Datasheet

GOT1 antibody (orb19657)

Description	Goat polyclonal to GOT1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GOT1 (aa 22-35)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 40kDa band observed in Human Liver lysates (calculated MW of 46.2kDa according to NP_002070.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Berk et al, Proc Natl Acad Sci U S A. 1990 May;87(9):3484-8; PMID: 2185471). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	TADFREDPDPKVN
Clonality	Polyclonal
MW	46.2
Entrez	2805
Dilution Range	ELISA: 1:16000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19657
(0.5ug/ml)
staining of
Human Li...