


Product Datasheet

HADH antibody (orb19669)

Description	Goat polyclonal antibody to HADH
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HADH / HADHSC
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Protein Sequence	YERGDASKEDID
Clonality	Polyclonal
MW	34.3
Entrez	3033
Dilution Range	ELISA: 1:64000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



orb19669
(0.02ug/ml)
staining of
Human K...



HEK293
overexpressing
HADH
(RC201752)
wi...