

Product Datasheet

YWHAQ antibody (orb19672)

Description	Goat polyclonal antibody to YWHAQ
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	14-3-3 tau / YWHAQ
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 28kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 27.8kDa according to NP_006817.1). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	DDRKQTIDNSQ
Clonality	Polyclonal
MW	27.8
Entrez	10971
Dilution Range	ELISA: 1:4000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb19672
(0.03ug/ml)
staining of
Human B...