


Product Datasheet

WISP1 antibody (orb19683)

Description	Goat polyclonal antibody to WISP1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	WISP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Protein Sequence	ESYPDFSEIAN
Clonality	Polyclonal
MW	40.3; 30.7; 13.5
Entrez	8840
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19683
(0.1ug/ml)
staining of
HEK293 c...



HEK293
overexpressing
WISP1
(RC214390)
W...