



Product Datasheet

PRRX1 antibody (orb19685)

Description	Goat polyclonal antibody to PRRX1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRRX1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 28kDa band observed in Human Heart lysates (calculated MW of 27.3kDa according to NP_073207.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	GSDTPQQDNDQ
Clonality	Polyclonal
MW	24.4; 27.3
Entrez	5396
Dilution Range	ELISA: 1:2000, WB: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ut...



Western blot analysis of LNCaP cell lysa...



Western blot analysis of Mouse Heart ysa...