


Product Datasheet

SOX15 antibody (orb19709)

Description	Goat polyclonal antibody to SOX15
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SOX15
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 24kDa band observed in Human Skeletal Muscle lysates (calculated MW of 25.3kDa according to NP_008873.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	CSLPQSDPRLQGE
Clonality	Polyclonal
MW	25.3
Entrez	6665
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19709
(0.1ug/ml)
staining of
Human Mu...