


Product Datasheet

ACVR1 antibody (orb19713)

Description	Goat polyclonal antibody to ACVR1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ACVR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Testis shows preferential staining of cytoplasm in primary spermatocytes.. Recommended concentration, 4-6ug/ml.WB: Approx 60kDa band observed in Human Umbilical Cord lysates (calculated MW of 57.2kDa according to NP_001096.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	RKFKRRNQERLNPRD
Clonality	Polyclonal
MW	57.2
Entrez	90
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Umbilical...



Immunohistochemical staining of Human Te...