

Product Datasheet

**Chicken anti Rabbit IgG (H+L), absorbed with
human and mouse IgG (orb1972356)**

Description

Chicken anti Rabbit IgG (H+L),

IgG

Species/Host

Gallus

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Form/Appearance

Phosphate buffered saline, pH 7.2
with 0.075% sodium azide

Concentration

See vial for concentration.

Storage

Storage: Product should be
stored at -20°C. Aliquot to avoid
freeze/thaw cycles

Note

For research use only

Source

Chickens were immunized with
highly purified rabbit IgG isolated
from normal pooled rabbit serum,
in Freund's adjuvant. Eggs were
collected, the IgY fraction
prepared and the anti rabbit IgG
immunoaffinity purified from a
solid phase immunoabsorbant.
Anti Rabbit IgG was further
absorbed against mouse and
human IgG.

**Hazard
Information**

This product is intended FOR
RESEARCH USE ONLY, and FOR
TESTS IN VITRO, not for use in
diagnostic or therapeutic
procedures involving humans or
animals.

Expiration Date

12 months from date of receipt.