

Product Datasheet

EPB41L4A antibody (orb19725)

Description

Goat polyclonal antibody to EPB41L4A

Species/Host

Goat

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

NBL4

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Preliminary experiments gave an approx 100kDa band in Human Brain (Cerebellum) lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 79.0kDa according to NP_071423.3. The 100kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

KNENSPDTQRSKSH

Clonality

Polyclonal

MW

79

Entrez

64097

Dilution Range

ELISA: 1:8000, WB: 1 µg/ml

Expiration Date

12 months from date of receipt.