

Product Datasheet

Anti-PPP1R37 Antibody (orb1972545)

Description

Biorbyt Anti-PPP1R37 Antibody catalog #

	applications. This anti...
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	E.coli-derived human PPP1R37 recombinant protein (Position: R66-R390). Human PPP1R37 shares 97.5% and 97.2% amino acid (aa) sequence identity with mouse and rat PPP1R37, respectively.
Form/Appearance	Lyophilized
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Clonality	Polyclonal
Uniprot ID	075864
Dilution Range	Western blot, 0.25-0.5 µg/ml, Human Flow Cytometry (Fixed), 1-3 µg /1x10 ⁶ cells, Human ELISA, 0.1-0.5 µg/ml, Human
Expiration Date	12 months from date of receipt.