

Product Datasheet

RPE65 antibody (orb19749)

Description	Goat polyclonal to RPE65
Species/Host	Goat
Reactivity	Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	RPE65
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60kDa band observed in Rat Retina lysates (calculated MW of 60.9kDa according to Human NP_000320.1 and 61.0kDa according to Rat NP_446014.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	EVKKNARKAPQPE
Clonality	Polyclonal
MW	60.9
Entrez	6121
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19749
(0.5ug/ml)
staining of
Rat Reti...