

Product Datasheet

PAM antibody (orb19751)

Description	Goat polyclonal to PAM
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PAM / PAL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise all reported isoforms (NP_000910.2; NP_620121.1; NP_620176.1 and NP_620177.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 100kDa band observed in Human Kidney lysates (calculated MW of 101kDa according to preproprotein NP_620121.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DQEKEDDGSESEE
Clonality	Polyclonal
MW	108; 101; 96.3; 98.8
Entrez	5066
Dilution Range	ELISA: 1:4000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19751
(1ug/ml)
staining of
Human Kidn...