

Product Datasheet

human CD29/PerCP Antibody (orb1975384)

Description	human CD29/PerCP Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD29
Preservatives	0.01M TBS(pH=7.4)
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
Uniprot ID	P05556
Expiration Date	12 months from date of receipt.