

Product Datasheet

PGC/Biotin Antibody (orb1975515)

Description	PGC/Biotin Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Recombinant human Pepsinogen II
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P20142
Expiration Date	12 months from date of receipt.