

Product Datasheet

TLR7/Biotin Antibody (orb1975616)

Description	TLR7/Biotin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant human TLR7 protein (24-348/1049aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Uniprot ID	Q9NYK1
Expiration Date	12 months from date of receipt.