

Product Datasheet

AIF1/Biotin Antibody (orb1975628)

Description	AIF1/Biotin Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Recombinant mouse AIF1 protein (2-147/147aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Uniprot ID	O70200
Expiration Date	12 months from date of receipt.