


Product Datasheet


CSNK1E antibody (orb19760)

Description	Goat polyclonal antibody to CSNK1E
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	CSNK1E
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: In transfected HEK293 transiently expressing CSNK1E a band of approx. 49kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	PASQTSVPFDHLGK
Clonality	Polyclonal
MW	47.3
Entrez	1454
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



HEK293
overexpressing
CSNK1E (RC202436)

...



Immunohistochemical
staining of Human
To...