

Product Datasheet

APOBEC3A antibody (orb19763)

Description

Goat polyclonal antibody to APOBEC3A

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Target

Phorbolin 1 / APOBEC3A

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:8000. WB: Preliminary experiments gave an approx 28kDa band in Human Spleen, Tonsil and Peripheral Blood Mononucleocyte lysates after 0.1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 23.0kDa according to NP_663745.1. The 28kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

TSNFNNGIGRHKTY

Clonality

Polyclonal

MW

23

Entrez


200315

Dilution Range


ELISA: 1:8000, WB: 0.1-1ug/ml, IHC-P: 3-6 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of Human To...



Western blot analysis of Human Spleen (A...



Flow cytometric analysis of HEK293 cells...