


Product Datasheet

Kcnj11 antibody (orb19764)

Description	Goat polyclonal antibody to Kcnj11
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Kcnj11 / Kir6.2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.IHC: In paraffin embedded Human Pancreas shows variable staining across the islet of Langerhans. Recommended concentration, 3-5ug/ml.WB: Approx 45kDa band observed in Human Skeletal Muscle lysates (calculated MW of 43.6kDa according to NP_034732.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	ERRARFVSKKGNC
Clonality	Polyclonal
MW	43.6
Entrez	3767
Dilution Range	ELISA: 1:64000, WB: 0.01-0.03 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



orb19764
(0.01ug/ml)
staining of
Human M...



orb19764
(3.8ug/ml)
staining of
paraffin...