


Product Datasheet

GOT2 antibody (orb19767)

Description	Goat polyclonal antibody to GOT2
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GOT2 (aa 295 to 306)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 45kDa band observed in Human Kidney and Liver lysates (calculated MW of 47.5kDa according to NP_002071.2). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	CKDADEAKRVES
Clonality	Polyclonal
MW	47.5
Entrez	2806
Dilution Range	ELISA: 1:32000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb19767
(0.03ug/ml)
staining of
Human K...