

Product Datasheet

41704 antibody (orb19771)

Description	Goat polyclonal antibody to 41704
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TEB4 / MARCH-VI (aa 701 to 712)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 100kDa band observed in Human Lung lysates (calculated MW of 103kDa according to NP_005876.2). Recommended concentration: 0.2-0.5 µg/ml. An additional band of unknown identity was also consistently observed at 18kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported
Protein Sequence	QGRRVIFQKVKE
Clonality	Polyclonal
MW	103
Entrez	10299
Dilution Range	ELISA: 1:16000, WB: 0.2-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lung lysa...