



Product Datasheet

CDKN2A Antibody (orb197810)


Description	Mouse monoclonal antibody to CDKN2A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Purified recombinant fragment of human CDKN2A (AA: 1-156) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	16.5kDa
Entrez	1029
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 (Lane 1)...



Western blot analysis of human CDKN2A us...



Western blot analysis of Hela (Lane 1), ...