



Product Datasheet

EIF2AK2 Antibody (orb197815)


Description	Mouse monoclonal antibody to EIF2AK2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Purified recombinant fragment of human EIF2AK2 (AA: 329-551) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	62kDa
Entrez	5610
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 (Lane 1), ...



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of A431 cells l...