



## Product Datasheet

### PAX3 Antibody (orb197831)


|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to PAX3                                               |
| <b>Species/Host</b>        | Mouse                                                                           |
| <b>Reactivity</b>          | Human, Mouse                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | ELISA, WB                                                                       |
| <b>Immunogen</b>           | Purified recombinant fragment of human PAX3 (AA: 142-203) expressed in E. Coli. |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                                    |
| <b>Concentration</b>       | 1mg/ml                                                                          |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | Mouse IgG1                                                                      |
| <b>Clonality</b>           | Monoclonal                                                                      |
| <b>MW</b>                  | 53kDa                                                                           |
| <b>Entrez</b>              | <b>5077</b>                                                                     |
| <b>Dilution Range</b>      | WB: 1/500 - 1/2000; ELISA: 1/10000                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human PAX3 usin...



Western blot analysis of HEK293 (Lane 1)...