


Product Datasheet

PNPLA2 antibody (orb19789)

Description	Goat polyclonal antibody to PNPLA2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PNPLA2 / ATGL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 55kDa band observed in Human Adipose lysates (calculated MW of 55.3kDa according to NP_065109.1). Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	PARPHGPEDKDQ
Clonality	Polyclonal
MW	55.3
Entrez	57104
Dilution Range	ELISA: 1:8000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



orb19789
(0.2ug/ml)
staining of
Human Ad...