


Product Datasheet

CES1 antibody (orb19790)

Description	Goat polyclonal antibody to CES1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CES1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: In paraffin embedded Human Lung shows staining of lymphocytes in the vasculature.. Recommended concentration, 2-3ug/ml.WB: Approx 55-60kDa band observed in Human Liver and Ileum lysates (calculated MW of 62.6kDa according to NP_001020366.1). In transfected HEK293 transiently expressing CES1 a band of approx. 65kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	NTQAAQKLKDKE
Clonality	Polyclonal
MW	62.6; 62.5; 62.4
Entrez	1066
Dilution Range	ELISA: 1:8000, WB: 0.03-0.1 µg/ml, IHC-P: 2-3µg/ml
Expiration Date	12 months from date of receipt.



orb19790
(0.03ug/ml)
staining of
Human L...



orb19790
(2.5ug/ml)
staining of
paraffin...