

Product Datasheet

Integrin beta 1 antibody (orb197998)

Description

Integrin beta 1 antibody

Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Gallus, Human, Mouse, Porcine, Rat, Sheep

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Integrin beta 1. Please contact us for the exact immunogen sequence. The peptide is available as orb374836.

Target

Integrin beta 1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

88 kDa

Uniprot ID

[P05556](#), [P09055](#), [P49134](#)

NCBI

[391988](#), [002202](#)

Entrez


[3688](#)

Dilution Range


IHC-P:1:100

Expiration Date


12 months from date of receipt.



IHC-P image of rat kidney tissue using I...



Immunohistochemical staining of pig hear...



Immunohistochemical staining of paraffin...