


Product Datasheet

ESRRG antibody (orb19800)

Description	Goat polyclonal antibody to ESRRG
Species/Host	Goat
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ESRRG
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx. 50kDa band observed in Human Brain (Cerebellum) and Breast lysates (calculated MW of 51.3kDa according to NP_NP_001429.2). In transfected HEK293 transiently expressing ESRRG a band of approx. 52kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	CRMSNKDRHIDS
Clonality	Polyclonal
MW	51.3; 48.6; 48.6
Entrez	2104
Dilution Range	ELISA: 1:16000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



orb19800
(0.2ug/ml)
staining of
Human Br...



HEK293
overexpressing
ESRRG
(RC218233)
a...