

Product Datasheet

EDG2 (orb198007)

Description	EDG2
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol
Form/Appearance	Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	For use as a positive control with Biorbyt's EDG2 (LPA1) polyclonal antibody (Cat. No. C174P) for Western blot analysis.
Uniprot ID	Q92633
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.