

Product Datasheet

EDG6 (orb198010)

Description	EDG6
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol
Form/Appearance	Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	For use as a positive control with Biorbyt's EDG6 (S1P4) polyclonal antibody (Cat. No. X1182P) for Western blot analysis.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.