

Product Datasheet

CD41 antibody (orb198064)

Description	Provided as solution in phosphate buffered saline with 0.08% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Flow Cytometry.
Isotype	IgG1
Clonality	Monoclonal
Source	96.2C1 was generated by immunizing mice with human PBL's.
Uniprot ID	P28906
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.