

Product Datasheet

IGFBP4 antibody (orb19808)

Description	Goat polyclonal antibody to IGFBP4
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IGFBP4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 30kDa band observed in Human Liver lysates (calculated MW of 27.9kDa according to NP_001543.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KGELDCHQLADSFRE
Clonality	Polyclonal
MW	27.9
Entrez	3487
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19808
(1ug/ml)
staining of
Human Live...