

Product Datasheet

HAND1 antibody (orb19812)

Description

HAND1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human HAND1. Please contact us for the exact immunogen sequence. The peptide is available as orb393983.

Target

HAND1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

23 kDa

Uniprot ID

[P10683](#)

NCBI

[015973](#), [057057](#)

Entrez



[51083](#)

Dilution Range


IHC-P: 1:100-500, WB: 1:200-800

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of mouse ovary...



WB analysis of rat ovary (lane 1), rat heart...



IHC-P image of pig small intestine tissue...